

# Pipette Controller

Levo



Compress the suction bellows before attaching the pipette

Button to expel residual liquid when using blow-out pipets

Filling and delivery are easily controlled with the small sensitive lever

The filter and adapter can be integral and replaceable

### Ordering information

Cat No.	Description
740100010000	Levo - Blue
740100020000	Levo - Red
740100030000	Levo - Green
740100040000	Levo - Yellow
740100050000	Levo - Pink

### Accessories

Cat No.	Description
17000110	Filter 3.0um, 1pc each pk

\* O.E.M. option also available.  
\* Design & Specification of Levo Pipette Controller & its accessories are subject to change without prior notice.

### Features:

- ❑ Comfortable and simple to use
- ❑ Precise pipetting control
- ❑ Robust and lightweight
- ❑ 1-100ml volumetric or serological pipettes can be used
- ❑ Uses integral and replaceable 3um filter
- ❑ Easy to maintain and clean
- ❑ Five colors
- ❑ Autoclavable silicone pipette holder

### Dragon Laboratory Instruments Limited

F3, Block 12, No.28 Yuhua Road, Shunyi District, Zone B,  
Tian Zhu Airport Industrial Park, Beijing, 101318 P. R. China  
Tel: 86-10-85653451 85653452 85653453  
Fax: 86-10-85653383  
E-mail: info@dragon-lab.com  
Service: service@dragon-lab.com  
Service Hotline: 400-1100-556  
Website: Http://www.dragon-lab.com

### US office

**DRAGONLAB, LLC**  
500 Four Rod Road, Suite 122, Berlin, CT 06037 USA  
Tel: 1-860-828-5289  
Fax: 1-860-828-5389  
E-mail: info@dragonlabusa.com  
Website: Http://www.dragon-lab.com